Analytix

Issue 3, Issue 4, Issue 5 • 2012, Issue 1 • 2013

Double Accreditation within Sigma-Aldrich®

SIGMA-ALDRICH®

Double Accreditation within Sigma-Aldrich®

Michael Weber
Manager Marketing and R&D
Analytical Standards and Reagents

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Kind regards,

Michael Weber
Manager Marketing and R&D
Analytical Standards and Reagents
michael.weber@sial.com

* With kind permission from Swiss Federal
Institute of Metrology (METAS)

Analytix is published five times per year by Sigma-Aldrich Chemie GmbH,
MarCom Europe, Industriestrasse 25, CH-9471 Buchs SG, Switzerland
Publisher: Sigma-Aldrich Marketing Communications Europe
Publication: Michael Jeitziner
Editor: Isabell Davies-Falch

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Alan Nichols, Manager R&D and Marketing alan.nichols@sial.com

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Figure 1 ISO/IEC 17025 and ISO Guide 34 Marks from the Swiss Accreditation Service (SAS) and the US ANSI-ASQ National Accreditation Board (ACLASS). Both SAS and ACLASS accreditation certificates are accepted worldwide.

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Table 1 Sigma-Aldrich has four sites which are doubly accredited and are producing certified reference materials.

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| Cat. no. | Brand | Description | Package size |
|----------|------------|---|--------------|
| CRM48392 | Supelco | Triazine Pesticide Standard, TraceCERT ®, 100 µg/mL each component in methanol Ametryn, atrazine, prometon, prometryn, propazine, simazine, terbutryn | 1 mL |
| CRM7505 | Supelco | HC BTEX/MTBE, TraceCERT , 2000 µg/mL each component in methanol Benzene, toluene, ethylbenzene, o-, m-, & p-xylene, methyl tert-butyl ether | 1 mL |
| CRM47177 | Supelco | Formaldehyde-DNPH, TraceCERT , 100 µg/mL in acetonitrile | 1 mL |
| CRM48902 | Supelco | EPA 8270 semi-volatile internal standard mix, TraceCERT , 2000 µg/mL each component 2 mL in methylene chloride Acenaphthene-d ₁₀ , chrysene-d ₁₂ , 1,4-dichlorobenzene-d ₄ , naphthalene-d ₈ , perylene-d ₁₂ , phenanthrene-d ₁₀ | 1 mL |
| CRM40071 | Supelco | Benzo(a)pyrene, TraceCERT , 1000 µg/mL each component in acetonitrile | 1 mL |
| 42438 | Fluka | Benzyl butyl phthalate, TraceCERT | 50 mg |
| 69311 | Fluka | Chinine, TraceCERT | 50 mg |
| 43540 | Fluka | Diisobutyl phthalate, TraceCERT | 100 mg |
| 51791 | Fluka | Acetate standard for IC, TraceCERT , 1000 mg/L acetate in water | 100 mL |
| 42071 | Fluka | Manganese standard for ICP, TraceCERT , 10 g/L manganese in 5% nitric acid | 100 mL |
| 83787 | Fluka | Rhenium standard for AAS, TraceCERT , 1 g/L rhenium in 2% nitric acid | 100 mL |
| PHR1238 | RTC | Ergocalciferol (vitamin D ₂) | 500 mg |
| PHR1185 | RTC | Amlodipine besylate | 1 g |
| PHR1196 | RTC | Nitrofurazone (nitrofur) | 1 g |
| PHR1235 | RTC | Retinyl palmitate (vitamin A palmitate) | 1 g |
| PHR1178 | RTC | Budesonide | 500 mg |
| B-050 | Cerilliant | Buphedrone ephedrine metabolite HCl, 1.0 mg/mL (as free base) in methanol | 1 mL |
| T-077 | Cerilliant | 3,3',5-Triiodo-L-thyronine- ¹³ C ₆ (T3- ¹³ C ₆), 100 µg/mL in 0.1N NH ₃ in methanol | 1 mL |
| A-096 | Cerilliant | Aldosterone, 100 µg/mL in acetonitrile | 1 mL |
| Z-009 | Cerilliant | Zopiclone N-oxide, 100 µg/mL in methanol | 1 mL |
| G-015 | Cerilliant | Ginseng ginsenosides mix, 100 µg/mL of each component in methanol | 1 mL |

Table 2 Recent certified reference materials released by Sigma-Aldrich

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| | |
|---|---|
| BAM – Federal Institute for Materials Research and Testing, Germany | MC-TIMS, double IDMS |
| CENAM – Centro Nacional de Metrología, Mexico | ICP-OES, one-point-calibration + IS |
| GUM – Central Office of Measures, Poland | ICP-OES, one-point-calibration |
| HKGL – Government Laboratory, Hong Kong | ICP-OES, one-point-calibration + IS |
| INM – National Institute of Metrology, Romania | Q-ICP-MS, calibration curve |
| INMETRO – National Institute of Metrology, Standardization and Industrial Quality, Brazil | ICP-OES, calibration curve |
| INTI – Instituto Nacional de Tecnología Industrial, Argentina | FAAS and ICP-OES, one-point-calibration |
| KRISS – Korea Research Institute of Standards and Science | ICP OES, one-point-calibration + IS |
| LGC – LGC Ltd., United Kingdom | MC-ICP-MS, double IDMS and Q-ICP-MS, one-point-calibration + IS |
| LNE – Laboratoire National de Métrologie et d'Essais, France | Q-ICP-MS, double IDMS and titrimetry |
| NIM – National Institute of Metrology, P. R. China | Q-ICP-MS, bracketing + IS |
| NIST – National Institute of Standards and Technology, USA | ICP-OES, one-point-calibration + IS |
| NMIA – National Measurement Institute, Australia | HR-ICP-MS, double IDMS |
| NMIJ – National Metrology Institute of Japan | titrimetry |
| NMISA – National Metrology Institute of South Africa | ICP-OES, bracketing + IS |
| PTB – Physikalisch-Technische Bundesanstalt, Germany | ICP-OES, bracketing + IS |
| SMU – Slovak Institute of Metrology | coulometric titrimetry and titrimetry |
| TUBITAK UME – TUBITAK National Metrology Institute, Turkey | Q-ICP-MS, four-point-calibration + IS and Q-ICP-MS, double IDMS |
| VNIIM – D. I. Mendeleev Institute for Metrology, Russian Federation | ICP-OES, one-point-calibration |

Table 1 Participants of CCQM-K87 in alphabetical order of their acronyms and technique used for measurements

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
Figure 1 Important steps to successful qNMR measurements

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
Figure 2 qNMR Spectrum of Acenaphthene and CRM Dimethylsulfone

| | |
|---------------------|--|
| P_{Sample} | Purity of sample as mass fraction (Acenaphthene, in % g/g) |
| P_{Ref} | Purity of reference as mass fraction (Dimethylsulfone, in % g/g) |
| I_{Sample} | Integral of the sample signal (Acenaphthene = 112.2648) |
| I_{Ref} | Integral of the reference signal (Dimethylsulfone = 100.0000) |
| N_{Sample} | Number of sample protons (Acenaphthene = 4) |
| N_{Ref} | Number of reference protons (Dimethylsulfone = 6) |
| M_{Sample} | Molecular mass of the sample (Acenaphthene = 154.21) |
| M_{Ref} | Molecular mass of the reference (Dimethylsulfone = 94.13) |
| m_{Sample} | Mass of sample (Acenaphthene) |
| m_{Ref} | Mass of reference (Maleic acid) |

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References

- [1] Analytix No. 3 (2010), "New generation of organic CRM", pp4-6.
- [2] F. Malz, H. Jancke, J. Pharm. Biomed. Anal. 38 (2005), pp813-23.
- [3] Eurachem/CITAC Guide, 2nd Ed. (2000), "Quantifying uncertainty in analytical measurement."
- [4] Weber M.; BERM-13, oral presentation on High Performance qNMR (2012).

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
Figure 1 Chemical structure of L-Methionine sulfoximine and methyl laurate

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
Figure 2 Chemical structures of the sulfonic acid esters

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Table 1 USP<232> analytes and permissible daily exposure (PDE) for pharmaceutical products as oral uptake. (n) = not a safety concern

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| Cat. no. | Product | Description |
|----------|--|--|
| 89395 | USP 232 Element Impurities, Standard 1 | 1.5 mg/L: As 5 mg/L: Pb 15 mg/L: Hg 25 mg/L: Cd 100 mg/L: Mo, V 500 mg/L: Ni 1000 mg/L: Cu in 10% in nitric acid, 100 mL FEP bottle |
| 39967 | USP 232 Element Impurities, Standard 2 | 100 mg/L: Ir, Os, Pd, Pt, Rh, Ru in 10% hydrochloric acid, 100 mL FEP bottle |
| 02650 | Nitric acid | ~65%, <i>TraceSELECT</i> ® Ultra, for ultra trace analysis, 250 mL and 1 L FEP bottles |
| 96208 | Hydrochloric acid | 30%, <i>TraceSELECT</i> Ultra, for ultra trace analysis, 250 mL and 1 L FEP bottles |
| 02658 | Hydrofluoric acid | 49%, <i>TraceSELECT</i> Ultra, for ultra trace analysis, 250 mL and 1 L FEP bottles |
| 77239 | Sulfuric acid | 95%, <i>TraceSELECT</i> Ultra, for ultra trace analysis, 250 mL and 1 L FEP bottles |
| 14211 | Water | <i>TraceSELECT</i> Ultra, ACS reagent, for ultra trace analysis, 1 L HDPE bottle |

Table 2 Range of Fluka® branded products that are used with USP <233>. Please find the complete list of inorganic trace analysis products at sigma-aldrich.com/traceanalysis

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| | Cat. No. | Brand | Product | Package Size |
|--|----------|--------|--|--------------|
| | 64397 | Fluka® | α-Amylcinnamaldehyde, mixture of cis and trans | 100 mg |
| | 50375 | Fluka | Bis(methylthio) methane | 100 mg |
| | 53877 | Fluka | 2,6-Dimethoxyphenol | 100 mg |
| | 76048 | Fluka | Dimethyl sulfide | 100 mg |
| | 44527 | Fluka | 4,5-Dimethyl-3-hydroxy-2,5-dihydrofuran-2-one | 100 mg |
| | 92255 | Fluka | δ-Dodecalactone | 100 mg |
| | 55821 | Fluka | Ethyl 3-hydroxybutyrate | 100 mg |
| | 94185 | Fluka | Ethyl maltol | 100 mg |
| | 90728 | Fluka | Ethylpyrazine | 100 mg |

Table 1 New Analytical Standards for the Analysis of Flavors and Fragrances

Dr. Matthias Nold, Product Manager, Analytical Standards matthias.nold@sial.com

| | Cat. No. | Brand | Product | Package Size |
|--|----------|-------|----------------------------------|--------------|
| | 75034 | Fluka | 2,3-Hexanedione | 100 mg |
| | 91316 | Fluka | cis-3-Hexen-1-ol | 100 mg |
| | 53407 | Fluka | 5-(Hydroxymethyl) furfural | 100 mg |
| | 79467 | Fluka | Linalyl formate | 100 mg |
| | 77938 | Fluka | Methyl 3-hydroxyhexanoate | 100 mg |
| | 80133 | Fluka | Methyl cyclopentenolone | 100 mg |
| | 69193 | Fluka | Methyl nicotinate | 100 mg |
| | 62535 | Fluka | Methyl p-tert-butylphenylacetate | 100 mg |
| | 61352 | Fluka | 4-Methyl-5-thiazolethanol | 100 mg |
| | 18073 | Fluka | 2-Methylanisole | 100 mg |
| | 30761 | Fluka | 2-Methylbutyl acetate | 100 mg |
| | 49597 | Fluka | 5-Methylquinoxaline | 100 mg |
| | 44638 | Fluka | 3-(Methylthio) propionaldehyde | 100 mg |
| | 75041 | Fluka | Neohesperidin dihydrochalcone | 100 mg |
| | 51573 | Fluka | γ-Octalactone | 100 mg |
| | 76046 | Fluka | Piperidine | 100 mg |
| | 94474 | Fluka | Sucrose octaacetate | 100 mg |
| | 95162 | Fluka | 2,3,5,6-Tetramethylpyrazine | 100 mg |
| | 76047 | Fluka | Tripropionin | 100 mg |

Reference

- [1] "Practical Analysis of Flavor and Fragrance Materials", Kevin Goodner, Russell Rouseff, John Wiley & Sons, 2011.

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Michael Jeitziner, Market Segment Manager Analytical Reagents & Standards michael.jeitziner@sial.com

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Table 1 Product information


Figure 1 The structure of E- and Z-Capsaicin


Figure 2 ^1H -NMR spectra of synthetically produced E-Capsaicin contaminated with Z-Capsaicin, Source: Syndeco GmbH, Switzerland

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
Figure 1 Molecular structure of deoxynivalenol (DON)

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
Figure 2 Separation of mycotoxins on Ascentis Express F5 (Positive Ion Mode). Red arrows indicate deoxynivalenol (DON) metabolites (the chromatogram was kindly provided by Enio Belotti, Water & Life Lab SRL, Entratico / Italy).

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
Figure 3 New mycotoxin standards brochure

| Cat. no. | Description | Concentration | Solvent | Package size |
|--|--|-----------------|------------------------------|--------------|
| Neats | | | | |
| 32927 | 3-acetyldeoxynivalenol | | | 5 mg |
| 32928 | 15-acetyldeoxynivalenol | | | 5 mg |
| 32943 | Deoxynivalenol | | | 5 mg |
| Solutions | | | | |
| 34132 | 3-acetyldeoxynivalenol | 100 µg/g | Acetonitrile | 2 mL |
| 34133 | 15-acetyldeoxynivalenol | 100 µg/g | Acetonitrile | 2 mL |
| 34129 | Acetyl-deoxynivalenol- D_3 | 100 µg/g | Acetonitrile | 2 mL |
| 46911 | Deoxynivalenol | 200 µg/g | Ethylacetate/Methanol (95/5) | 1 mL |
| 34155 | Deoxynivalenol- D_1 | 100 µg/g | Acetonitrile | 2 mL |
| 34135 | Deepoxy-deoxynivalenol | 50 µg/g | Acetonitrile | 2 mL |
| 34124 | Deoxynivalenol | 100 µg/g | Acetonitrile | 2 mL |
| 32911 | Deoxynivalenol-3-glucoside | 50 µg/g | Acetonitrile | 1 mL |
| ^{13}C-Isotopical Labeled | | | | |
| 32962 | 3-acetyldeoxynivalenol- $^{13}\text{C}_{17}$ | 25 µg/mL | Acetonitrile | 1 mL |
| 34128 | Deoxynivalenol- $^{13}\text{C}_{15}$ | 25 µg/mL | Acetonitrile | 1 mL |
| Mixture Solutions | | | | |
| 32926 | Trichothecene Mix (3-AcDON, DON, NIV, FusX, HT-2, T-2, DAS, ZON) | 10 µg/mL (each) | Acetonitrile | 1 mL |
| 34134 | B-Trichothecene Mix (DON, NIV, 3-AcDON, 15-AcDON) | 100 µg/mL | Acetonitrile | 2 mL |
| Reference Materials | | | | |
| BCR377 | Maize flour (Deoxynivalenol, blank) | <50 µg/kg | | 150 g |
| IRMM315 | 4-Deoxynivalenol | 25.1 µg/g | Acetonitrile | 4 mL |

Table 1 Analytical standards for deoxynivalenol (DON) metabolite analysis

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
Figure 1 ESI LC-MS spectra of lithium formate solution (Fluka®) in positive (upper) and negative ion mode (bottom). The regular pattern of the $[Li_{x+1}(HCOO)_x]^+$ and $[Li_x(HCOO)_{x+1}]^-$ cluster can be used to calibrate a mass range up to 1200 m/z , which covers most of the small molecules.

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Reference

[1] J. Emmert, Labor & More, 5 (2008), 56-57.

| Cat. No. | Brand | Description |
|------------|--------|---|
| 01886 | Fluka® | Lithium Formate Solution, Suitable for LC-MS, 10 mM LiOH in isopropanol/water 1:1 (+0.2% HCOOH) |
| 97574 | Fluka | Sodium Formate Solution, Suitable for LC-MS, 10 mM NaOH in isopropanol/water 1:1 (+0.2% HCOOH) |
| 43530 | Fluka | Reserpine Standard for LC-MS, analytical standard, for LC-MS |
| 56396 | Fluka | Caffeine, certified reference material, TraceCERT® |
| MSCAL1-1KT | Sigma | ProteoMass Peptide and Protein MALDI-MS Calibration Kit |
| D6003 | Fluka | Digoxin, analytical standard |


Figure 2 Influence of an aged column on the performance of a separation. 5 µl of a 5 ppm reserpine standard is injected and monitored over a column lifetime. Peak height and symmetry decreases. UHPLC separations especially react sensitively to instrumental errors (extra volume, valve malfunction, clogging, or performance loss in the stationary phase).


Figure 3 Acquired and calculated mass spectrum of reserpine. The difference of the peak position is only 0.7 ppm. The resolution is 11000 in standard resolution mode. This agrees with the specification of the manufacturer and indicates an optimum setup.

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Matthias Nold, Product Manager, Analytical Standards matthias.nold@sial.com

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(continued on page 7)


Figure 1 Molecular structures of ^{13}C ring labeled parabens

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
Figure 2 Chromatogram of parabens and ^{13}C -labeled internal standards

GC Conditions:

instrument: Agilent GC 7890A
 detector: Agilent MS 7000 QQQ
 column: HP-5MS, 5% phenyl polymethylsiloxane, 25 m x 0.25 mm, 0.25 μm (L x ID, FT)
 oven program: 760 $^{\circ}\text{C}$ (1 min), 10 $^{\circ}\text{C}/\text{min}$ to 200 $^{\circ}\text{C}$, 15 $^{\circ}\text{C}/\text{min}$ to 280 (10 min)
 carrier gas: helium, 1.2 mL/min
 injection: 2.0 μL splitless
 injector: Gerstel CIS
 injector temperature program: 60 $^{\circ}\text{C}$, 12 $^{\circ}\text{C}/\text{to}$ 280 $^{\circ}\text{C}$ (2 min)

(continued on page 8)

| Compound | Cat. no. | Concentration | Package size | Cat. no. | Concentration | Package size |
|---|----------|---------------------|--------------|----------|---------------------|--------------|
| Methylparaben (Ring- ¹³ C ₆) | 32557 | 50 µg/mL in acetone | 1 mL | 32556 | 10 µg/mL in acetone | 5 mL |
| Ethylparaben (Ring- ¹³ C ₆) | 32559 | 50 µg/mL in acetone | 1 mL | 32558 | 10 µg/mL in acetone | 5 mL |
| n-Propylparaben (Ring- ¹³ C ₆) | 32562 | 50 µg/mL in acetone | 1 mL | 32561 | 10 µg/mL in acetone | 5 mL |
| n-Butylparaben (Ring- ¹³ C ₆) | 32564 | 50 µg/mL in acetone | 1 mL | 32563 | 10 µg/mL in acetone | 5 mL |

Table 1 ¹³C labeled paraben standard solutions

| Brand | Cat. no. | Description | Package size |
|--------|----------|---|--------------|
| Fluka® | PHR1012 | Methylparaben pharmaceutical secondary standard | 1 g |
| Fluka | PHR1011 | Ethylparaben pharmaceutical secondary standard | 1 g |
| Fluka | PHR1010 | n-Propylparaben pharmaceutical secondary standard | 1 g |
| Fluka | PHR1022 | n-Butylparaben pharmaceutical secondary standard | 1 g |

Table 2 Non labeled paraben standards

| Description | Components | Cat. no. | Concentration | Package size | Cat. no. | Concentration | Package size |
|--|--|----------|------------------------------------|--------------|----------|-----------------------------------|--------------|
| Paraben Internal Standard Mix-Solution | Methylparaben Ring- ¹³ C ₆ Ethylparaben Ring- ¹³ C ₆ n-Propylparaben Ring- ¹³ C ₆ n-Butylparaben Ring- ¹³ C ₆ | 32124 | 10 µg/mL each component in acetone | 5mL | 32125 | 50 µg/mL each component in Aceton | 1mL |
| Paraben Target Analyte Mix-Solution | Methylparaben Ethylparaben n-Propylparaben iso-Propylparaben n-Butylparaben iso-Butylparaben Benzylparaben | 32126 | 10 µg/mL each component in acetone | 5mL | 32127 | 50 µg/mL each component in Aceton | 1mL |

Table 3 Paraben Mix-Solutions

References

- [1] Routledge, E.J., Parker, J., Odum, J., Ashby, J., Sumpter, J.P. Some alkyl hydroxy benzoate preservatives (parabens) are estrogenic, *Toxicol. Appl. Pharmacol.* 153 (1998) 12–19.
- [2] Canosa, P., Rodriguez, I., Rubi, E., Cela, R. Determination of parabens and triclosan in indoor dust using matrix solid-phase dispersion and gas chromatography with tandem mass spectrometry, *Anal. Chem.* 79 (2007) 1675–1681.
- [3] Blanco, E., Maria del, C.C., Maria del, C.M., Cela, R. Combination of off-line solid-phase extraction and on-column sample stacking for sensitive determination of parabens and p-hydroxybenzoic acid in waters by non-aqueous capillary electrophoresis, *Anal. Chim. Acta* 647 (2009) 104–111.
- [4] Canosa, P., Rodriguez, I., Rubi, E., Bollain, M.H., Cela, R. Optimisation of a solid-phase microextraction method for the determination of parabens in water samples at the low ng per litre level, *J. Chromatogr. A* 1124 (2006) 3–10.
- [5] Ye, X., Bishop, A.M., Needham, L.L., Calafat, A.M. Automated on-line column switching HPLC–MS/MS method with peak focusing for measuring parabens, triclosan, and other environmental phenols in human milk, *Anal. Chim. Acta* 622 (2008) 150–156.
- [6] Janjua, N.R., Mortensen, G.K., Andersson, A.M., Kongshoj, B., Skakkebaek, N.E., Wulf, H.C. Systemic uptake of diethyl phthalate, dibutyl phthalate, and butyl paraben following whole-body topical application and reproductive and thyroid hormone levels in humans, *Environ. Sci. Technol.* 41 (2007) 5564–5570.
- [7] Ye, X., Bishop, A.M., Reidy, J.A., Needham, L.L., Calafat, A.M. Parabens as urinary biomarkers of exposure in humans, *Environ. Health Perspect.* 114 (2006) 1843–1846.
- [8] Darbre, P.D., Aljarrah, A., Miller, W.R., Coldham, N.G., Sauer, M.J. Concentrations of parabens in human breast tumours, *J. Appl. Toxicol.* 24 (2004) 5–13.
- [9] Kang, K.S., Che, J.H., et al Decreased sperm number and motile activity on the F1 offspring maternally exposed to butyl p-hydroxybenzoic acid (butyl paraben), *The Journal of Veterinary Medical Science* 64 (2002), 227–235.
- [10] Nieto, A., Borull, F., Marce, R.M., Pocurull, E. Determination of personal care products in sewage sludge by pressurized liquid extraction and ultra high performance liquid chromatography – tandem mass spectrometry, *J. Chromatogr. A* 1216 (2009) 5619–5625.
- [11] SCF 1996, Scientific Committee on Food, Opinion on p-hydroxybenzoic acid alkyl esters and their sodium salts, 25 February 1994, 35th series of reports of the SCF.
- [12] SCCS 2011, Scientific Committee on Consumer Safety, Opinion on Parabens, 22 March 2011, SCCS/1348/10.
- [13] SCF 2003, SCF/CS/ADD/CONS/53 Final: Statement of the Scientific Committee on Food on Parabens.

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Rudolf Köhling, Senior Scientist, R&D LC-MS Applications rudolf.koehling@sial.com

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
Figure 1 a) UV chromatograms of 5 injections of vitamin D₃ at 254 nm using 2 different grades of acetonitrile. b) UV chromatogram of 5 injections of 25-hydroxyvitamin D₃ repeating the experiment of vitamin D₃. The UV data indicates the robustness of both the chromatographic separation and the precision of the UHPLC-MS system.

| | |
|---|--|
| Mass spectrometer (Bruker microTOF-Q II) | tune_low.m, 45 psi (sprayer), 7 L/min (dry gas), 180 °C (dry gas temp), ESI(+) |
| HPLC (Dionex Ultimate 3000 RSLC) Column: | Supelco® Ascentis® Express C18, 2.1x50 mm or comparable column |
| Flow: | 0.6 mL/min |
| Gradient: | 0 min (50% C), 5 min (10% C) Post time: 2.5 min |
| Temperature: | 35 °C |
| Injection Volume: | 2 µL (vitamin D ₃)/4 µL (25-hydroxyvitamin D ₃) |
| UV Wave Length: | 254 nm |

Table 1 LC-MS Method for Vitamin D₃/Hydroxyvitamin D₃

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
Figure 2 a) EICs ($[M(\text{vitamin D}_3)+H]^+=385.4$ Da) of LC-MS Ultra acetonitrile and competitor acetonitrile. The loss of sensitivity concerning the isomer of vitamin D₃ at 3.68 min is remarkable. With LC-MS Ultra acetonitrile, this isomer is detectable and it is below the LOD with the lower quality solvent. b) EICs ($[M(25\text{-hydroxyvitamin D}_3)+H]^+=407.4$ Da) of repeated experiments using different quality grades of acetonitrile show the same ion suppression effect with the lower quality solvent.

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
Figure 3 UV chromatograms of different UHPLC/UHPLC-MS acetonitrile manufacturers tested with the UHPLC gradient method against LC-MS Ultra water. Even at 210 nm, there are only minor differences.


Figure 4 The base peak chromatograms (BPC) of the corresponding UV data in Figure 2 reveal the real differences and show a higher baseline and/or additional peaks compared to the reference solvents (LC-MS Ultra). Quality M and R are comparable although one is tested for LC-MS and the other is not.

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(continued on page 16)



Figure 5 Background of vitamin D₃/25-hydroxyvitamin D₃ chromatogram with low-quality acetonitrile. The spectrum is typical for alkali formate clusters. The high abundance of sodium can be responsible for this considerable ion suppression effect.

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
References

- [1] Pucci, V., Di Palma, S., Alfieri, A., Bonelli, F., Monteagudo, E. A Novel Strategy for Reducing Phospholipids-based Matrix Effect in LC-ESI-MS Bioanalysis by Means of HybridSPE. *J. Pharm. Biomed. Anal.* 2009, 50(5), 867–871.
- [2] Ismael, O. A., Zhang, T., Jenkins, R. G.; Karnes, H. T. Investigation of Endogenous Blood Plasma Phospholipids, Cholesterol and Glycerides that Contribute to Matrix Effects in Bioanalysis by Liquid Chromatography/Mass Spectrometry. *J. Chromatogr., B*, 2010, 878, 3303–3316.
- [3] Mano, N., Nozawa, M.; Sato, M.; Mori, M., Yamaguchi, H.; Kanda, K., Nogami, M.; Goto, J., Shimada, M. Identification and Elimination of Ion Suppression in the Quantitative Analysis of Sirolimus in Human Blood by LC/ESI-MS/MS. *J. Chromatogr., B*, 2011, 879, 968–974.
- [4] Annesley, T. Ion Suppression in Mass Spectrometry. *Clinical Chemistry*, 2003, 49, 7, 1041–1044.

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Shyam Verma, Market Segment Manager shyam.verma@sial.com

Contributed article

Melissa P Grella, Jessica Hoover, and Mark Shapiro
PharmaCore, Inc, High Point, NC 27265, USA

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
Figure 1 DMA comparison: headspace grade vs synthetic grade

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

**Figure 2** Example Method 1**Figure 3** Example Method 2

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References

- [1] United States Pharmacopoeia (USP), 31st Edition (2008), <467> Residual Solvents.
- [2] ICH Guidelines for Industry, Q3C Impurities: Residual Solvents, US Dept. of Health and Human Services Food and Drug Administration, Center for Drug Evaluation and Research (CDER), Center for Biologics Evaluation and Research (CBER), ICH December, 1997.

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
Product Table Ionic liquid matrices (ILMs) for MALDI-MS. Please find the complete product list at sigma-aldrich.com/maldi

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Product Table HYDRANAL® Standard Substances for KF Titration

| Cat. No. | Brand | Description |
|-------------------------|--------|--|
| Liquid Standards | | |
| 34849 | Fluka® | HYDRANAL-Water Standard 10.0 Liquid standard for volumetric KF titer determination. Water content 10.0 mg/g (1.0%), exact value stated on CoA. Tested against NIST SRM 2890 and traceable to SI-Unit (kg). |
| 34828 | Fluka | HYDRANAL-Water Standard 1.0 Liquid standard for coulometric KF titration. Water content 1.0 mg/g (0.1%), exact value stated on CoA. Tested against NIST SRM 2890. Used for validation and check of system and recovery rates. |
| 34847 | Fluka | HYDRANAL-Water Standard 0.1 Liquid standard for coulometric KF titration. Water content 0.1 mg/g (0.01%), exact value stated on CoA. Tested against NIST SRM 2890. Used for validation and check of system and recovery rates. |
| 34813 | Fluka | HYDRANAL-Standard 5.0 Liquid standard for volumetric KF titration. Water content 5.00 ± 0.02 mg/mL (5.91 ± 0.02 mg/g). |
| 34694 | Fluka | HYDRANAL-Water Standard Oil Mineral oil-based matrix standard, for direct coulometric and volumetric KF titration and indirect titration with KF oven. Water content in the low ppm range, exact value stated on CoA. Tested against NIST SRM 2890. |
| Solid Standards | | |
| 34696 | Fluka | HYDRANAL-Standard sodium tartrate dihydrate Primary standard reference material for volumetric KF titration. Water content ~15.66%, exact value stated on CoA. Tested against NIST SRM 2890 and traceable to SI-Unit (kg). This standard remains stable under normal conditions and does not lose or absorb moisture. This product is a finely powdered solid that dissolves relatively quickly; although its solubility in methanol is limited. |
| 34693 | Fluka | HYDRANAL-Water Standard KF-Oven 140–160 °C Solid standard for the control of KF oven using the coulometric or volumetric technique. Water content ~5%, exact value stated on CoA. Tested against NIST SRM 2890. This standard can be used to verify accuracy of results and tightness of a KF oven system. The function of the connected titration cell should be tested using a liquid HYDRANAL-Water Standard. |
| 34748 | Fluka | HYDRANAL-Water Standard KF-Oven 220–230 °C Solid standard for the control of KF oven, using the coulometric or volumetric technique. Water content 5.55 ± 0.05%, exact value stated on CoA. Tested against NIST SRM 2890. This standard can be used to verify accuracy of results and tightness of a KF oven system. The function of the connected titration cell should be tested using a liquid HYDRANAL-Water Standard. |

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©2013 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. FLUKA is a trademark of Sigma-Aldrich GmbH, registered in the US and other countries. SAFC brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at sigmaaldrich.com. Ascentis, Cerilliant, HYDRANAL, *TraceCERT* and *TraceSELECT* are registered trademarks of Sigma-Aldrich Co. LLC.

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